

— ABLYBIO, Help Your Research



PLCB3 Antibody

货号: **AYP4643**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	131kDa/138kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	OVCAR8,Mouse lung,Rat lung
细胞定位	Membrane,Peripheral membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human PLCB3.
------	---

靶点信息

研究背景	This gene encodes a member of the phosphoinositide phospholipase C beta enzyme family that catalyze the production of the secondary messengers diacylglycerol and inositol 1,4,5-triphosphate from phosphatidylinositol in G-protein-linked receptor-mediated signal transduction. Alternative splicing results in multiple transcript variants.
基因ID	5331
基因名	PLCB3
Swiss	Q01970 (https://www.uniprot.org/uniprotkb/Q01970/entry)
别名	PLCB3,PLCB3 Antibody,Phosphoinositide phospholipase C-beta-3,Phospholipase C-beta-3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)