

— ABLYBIO, Help Your Research



# STAT5B Antibody

货号: **AYP4647**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	90kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K562
细胞定位	cytosol,nucleoplasm
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human STAT5B.
------	--

## 靶点信息

研究背景	Carries out a dual signal transduction and activation of transcription. Mediates cellular responses to the cytokine KITLG/SCF and other growth factors. Mediates cellular responses to ERBB4. May mediate cellular responses to activated FGFR1, FGFR2, FGFR3 and FGFR4. Binds to the GAS element and activates PRL-induced transcription. Regulates the expression of milk proteins during lactation.
基因ID	6776,6777
基因名	STAT5A,STAT5B
Swiss	P42229P51692 ( <a href="https://www.uniprot.org/uniprotkb/P42229P51692/entry">https://www.uniprot.org/uniprotkb/P42229P51692/entry</a> )
别名	MGF,STAT5,stat5a/stat5b,STAT5B,STAT5B Antibody,STAT5A

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)