

— ABLYBIO, Help Your Research



# Fyn Antibody

货号: **AYP4658**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	54kDa/60kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431
细胞定位	Cell membrane,Cytoplasm,Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human Fyn.
------	---

## 靶点信息

研究背景	This gene is a member of the protein-tyrosine kinase oncogene family. It encodes a membrane-associated tyrosine kinase that has been implicated in the control of cell growth. The protein associates with the p 85 subunit of phosphatidylinositol 3-kinase and interacts with the fyn-binding protein. Alternatively spliced transcript variants encoding distinct isoforms exist.
基因ID	2534
基因名	FYN
Swiss	P06241 ( <a href="https://www.uniprot.org/uniprotkb/P06241/entry">https://www.uniprot.org/uniprotkb/P06241/entry</a> )
别名	FYN,SLK,SYN,p59-FYN,Fyn Antibody,Proto-oncogene Syn,Proto-oncogene c-Fyn,Src-like kinase

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)