

Integrin β 1 Antibody

货号: **AYP4666**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	87kDa/88kDa/91kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,U-87MG
细胞定位	Cell junction,Cell membrane,Cell projection,Cleavage furrow,Melanosome,Recycling endosome,Single-pass type I membrane protein,invadopodium membrane,lamellipodium,ruffle membrane,ruffle,sarcolemma
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Integrin β 1.
------	--

靶点信息

研究背景	Integrins are heterodimeric proteins made up of alpha and beta subunits. At least 18 alpha and 8 beta subunits have been described in mammals. Integrin family members are membrane receptors involved in cell adhesion and recognition in a variety of processes including embryogenesis, hemostasis, tissue repair, immune response and metastatic diffusion of tumor cells. This gene encodes a beta subunit. Multiple alternatively spliced transcript variants which encode different protein isoforms have been found for this gene.
基因ID	3688
基因名	ITGB1

Swiss	P05556
别名	CD29;FNRB;GPIIA;MDF2;MSK12;VLA-BETA;VLAB;Integrin beta 1;ITGB1;Integrin β 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn