

— ABLYBIO, Help Your Research



ACK1 Antibody

货号: **AYP4683**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	60kDa/114kDa/119kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,HT-29,MCF7,Raji
细胞定位	Cell junction,Cell membrane,Cytoplasmic side,Cytoplasmic vesicle,Cytoplasmic vesicle membrane,Endosome,Membrane,Nucleus,Peripheral membrane protein,adherens junction,clathrin-coated pit,clathrin-coated vesicle
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ACK1.
------	--

靶点信息

研究背景	This gene encodes a tyrosine kinase that binds Cdc42Hs in its GTP-bound form and inhibits both the intrinsic and GTPase-activating protein (GAP)-stimulated GTPase activity of Cdc42Hs. This binding is mediated by a unique sequence of 47 amino acids C-terminal to an SH3 domain. The protein may be involved in a regulatory mechanism that sustains the GTP-bound active form of Cdc42Hs and which is directly linked to a tyrosine phosphorylation signal transduction pathway. Several alternatively spliced transcript variants have been identified from this gene, but the full-length nature of only two transcript variants has been determined.
基因ID	10188
基因名	TNK2
Swiss	Q07912 (https://www.uniprot.org/uniprotkb/Q07912/entry)
别名	TNK2,ACK,ACK-1,ACK1,p21cdc42Hs,ACK1 Antibody,Tyrosine kinase non-receptor protein 2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)