

— ABLYBIO, Help Your Research



# APC Antibody

货号: **AYP4700**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	300kDa/311kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,Mouse brain,Rat brain
细胞定位	Cell junction,Cell membrane,Cell projection,Cytoplasm,adherens junction,cytoskeleton,lamellipodium,ruffle membrane
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human APC.
------	---

## 靶点信息

研究背景	This gene encodes a tumor suppressor protein that acts as an antagonist of the Wnt signaling pathway. It is also involved in other processes including cell migration and adhesion, transcriptional activation, and apoptosis. Defects in this gene cause familial adenomatous polyposis (FAP), an autosomal dominant pre-malignant disease that usually progresses to malignancy. Disease-associated mutations tend to be clustered in a small region designated the mutation cluster region (MCR) and result in a truncated protein product.
基因ID	324
基因名	APC
Swiss	P25054 ( <a href="https://www.uniprot.org/uniprotkb/P25054/entry">https://www.uniprot.org/uniprotkb/P25054/entry</a> )
别名	APC,BTPS2,DP2,DP2.5,DP3,GS,PPP1R46,APC Antibody,Deleted in polyposis 2.5

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)