

— ABLYBIO, Help Your Research



# CTNND1 Antibody

货号: **AYP4705**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	65-72kDa/90-108kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,NIH/3T3,C6
细胞定位	Cell membrane,Cytoplasm,Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human CTNND1.
------	--

## 靶点信息

研究背景	This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length natures of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene.
基因ID	1500
基因名	CTNND1
Swiss	O60716 ( <a href="https://www.uniprot.org/uniprotkb/O60716/entry">https://www.uniprot.org/uniprotkb/O60716/entry</a> )
别名	CAS,CTNND,P120CAS,P120CTN,p120,p120(CAS),p120(CTN),CTNND1,CTNND1 Antibody,Cadherin-associated Src substrate,p120 catenin,KIAA0384

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)