

— ABLYBIO, Help Your Research



NUR77 Antibody

货号: **AYP4733**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200
理论分子量	33kDa/64kDa/65kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human NUR77.
------	---

靶点信息

研究背景	This gene encodes a member of the steroid-thyroid hormone-retinoid receptor superfamily. Expression is induced by phytohemagglutinin in human lymphocytes and by serum stimulation of arrested fibroblasts. The encoded protein acts as a nuclear transcription factor. Translocation of the protein from the nucleus to mitochondria induces apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	3164
基因名	NR4A1
Swiss	P22736 (https://www.uniprot.org/uniprotkb/P22736/entry)
别名	NR4A1,GFRP1,HMR,N10,NAK-1,NGFIB,NP10,NUR77,TR3,NUR77 Antibody,Early response protein NAK1,Nuclear hormone receptor NUR/77,Orphan nuclear receptor HMR,Orphan nuclear receptor TR3,ST-59,Testicular receptor 3,NAK1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)