

— ABLYBIO, Help Your Research



GEF H1 Antibody

货号: **AYP4737**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	108kDa/111kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C6
细胞定位	Cell junction,Cell projection,Cytoplasm,Cytoplasmic vesicle,Golgi apparatus,cytoskeleton,ruffle membrane,spindle,tight junction
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human GEF H1.
------	--

靶点信息

研究背景	Rho GTPases play a fundamental role in numerous cellular processes that are initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form complex with G proteins and stimulate rho-dependent signals. Alternatively spliced transcript variants encoding different isoforms have been identified.[provided by RefSeq, Jun 2009]
基因ID	9181
基因名	ARHGEF2
Swiss	Q92974 (https://www.uniprot.org/uniprotkb/Q92974/entry)
别名	GEF,GEF-H1,GEFH1,LFP40,P40,ARHGEF2,NEDMHM,GEF H1 Antibody,Guanine nucleotide exchange factor H1,Microtubule-regulated Rho-GEF,Proliferating cell nucleolar antigen p40,KIAA0651,GEF H1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)