

— ABLYBIO, Help Your Research



# FGR Antibody

货号: **AYP4804**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	59kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,U-937,Raji
细胞定位	Cell membrane,Cell projection,Cytoplasm,Cytoplasmic side,Lipid-anchor,Mitochondrion inner membrane, Mitochondrion intermembrane space,Peripheral membrane protein,cytoskeleton,cytosol,ruffle membrane
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human FGR.
------	---

## 靶点信息

研究背景	This gene is a member of the Src family of protein tyrosine kinases (PTKs). The encoded protein contains N-terminal sites for myristylation and palmylation, a PTK domain, and SH2 and SH3 domains which are involved in mediating protein-protein interactions with phosphotyrosine-containing and proline-rich motifs, respectively. The protein localizes to plasma membrane ruffles, and functions as a negative regulator of cell migration and adhesion triggered by the beta-2 integrin signal transduction pathway. Infection with Epstein-Barr virus results in the overexpression of this gene. Multiple alternatively spliced variants, encoding the same protein, have been identified.
基因ID	2268
基因名	FGR
Swiss	P09769 ( <a href="https://www.uniprot.org/uniprotkb/P09769/entry">https://www.uniprot.org/uniprotkb/P09769/entry</a> )
别名	FGR, SRC2, c-fgr, c-src2, p55-Fgr, p55c-fgr, p58-Fgr, p58c-fgr, FGR Antibody, Gardner-Rasheed feline sarcoma viral (v-fgr) oncogene homolog, Proto-oncogene c-Fgr

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)