

— ABLYBIO, Help Your Research



# ICK Antibody

货号: **AYP4805**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	34kDa/71kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SW620,Mouse liver,Mouse lung,Mouse brain
细胞定位	Cell projection,Cytoplasm,Nucleus,cilium,cilium basal body,cytoskeleton,cytosol
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human ICK.
------	---

## 靶点信息

研究背景	Eukaryotic protein kinases are enzymes that belong to a very extensive family of proteins which share a conserved catalytic core common with both serine/threonine and tyrosine protein kinases. This gene encodes an intestinal serine/threonine kinase harboring a dual phosphorylation site found in mitogen-activating protein (MAP) kinases. The protein localizes to the intestinal crypt region and is thought to be important in intestinal epithelial cell proliferation and differentiation. Alternative splicing has been observed at this locus and two variants, encoding the same isoform, have been identified.
基因ID	22858
基因名	ICK
Swiss	Q9UPZ9 ( <a href="https://www.uniprot.org/uniprotkb/Q9UPZ9/entry">https://www.uniprot.org/uniprotkb/Q9UPZ9/entry</a> )
别名	ICK,ECO,LCK2,MRK,ICK Antibody,Ciliogenesis associated kinase 1,Intestinal cell kinase,Laryngeal cancer kinase 2,MAK-related kinase,KIAA0936

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)