

— ABLYBIO, Help Your Research



MAK Antibody

货号: **AYP4813**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Rat lung
细胞定位	axoneme,centrosome,cytoplasm,mitotic spindle,nucleus,photoreceptor inner segment,photoreceptor outer segment
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human MAK.
------	---------------------------------------------

靶点信息

研究背景	The product of this gene is a serine/threonine protein kinase related to kinases involved in cell cycle regulation. Studies of the mouse and rat homologs have localized the kinase to the chromosomes during meiosis in spermatogenesis, specifically to the synaptonemal complex that exists while homologous chromosomes are paired. Mutations in this gene have been associated with ciliary defects resulting in retinitis pigmentosa 62. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]
基因ID	4117
基因名	MAK
Swiss	P20794 (https://www.uniprot.org/uniprotkb/P20794/entry)
别名	RP62,MAK,MAK Antibody,Male germ cell-associated kinase

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)