

— ABLYBIO, Help Your Research



# MSK2 Antibody

货号: **AYP4822**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	84kDa/85kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human MSK2.
------	--

## 靶点信息

研究背景	This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 non-identical kinase catalytic domains and phosphorylates various substrates, including CREB1 and ATF1. The encoded protein can also phosphorylate histone H3 to regulate certain inflammatory genes. Several transcript variants encoding different isoforms have been found for this gene.
基因ID	8986
基因名	RPS6KA4
Swiss	O75676 ( <a href="https://www.uniprot.org/uniprotkb/O75676/entry">https://www.uniprot.org/uniprotkb/O75676/entry</a> )
别名	RPS6KA4,MSK2,RSK-B,S6K-alpha-4,MSK2 Antibody,90 kDa ribosomal protein S6 kinase 4,Nuclear mitogen- and stress-activated protein kinase 2,Ribosomal protein kinase B

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))