

— ABLYBIO, Help Your Research



MYT1 Antibody

货号: **AYP4823**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	54kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse testis,Mouse thymus,Rat testis
细胞定位	cytosol,endoplasmic reticulum,Golgi apparatus,nucleolus,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human MYT1.
------	--

靶点信息

研究背景	This gene encodes a member of the serine/threonine protein kinase family. The encoded protein is a membrane-associated kinase that negatively regulates the G2/M transition of the cell cycle by phosphorylating and inactivating cyclin-dependent kinase 1. The activity of the encoded protein is regulated by polo-like kinase 1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	9088
基因名	PKMYT1
Swiss	Q99640 (https://www.uniprot.org/uniprotkb/Q99640/entry)
别名	MYT1,PPP1R126,MYT1 Antibody,PKMYT1,Myt1 kinase

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)