

— ABLYBIO, Help Your Research



# p70 S6 Kinase beta Antibody

货号: **AYP4824**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	16kDa/53kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SH-SY5Y,HL-60,SW620,MCF7,Mouse liver,Mouse kidney,Mouse spleen,Rat liver
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human p70 S6 Kinase beta.
------	--

## 靶点信息

研究背景	This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains a kinase catalytic domain and phosphorylates the S6 ribosomal protein and eukaryotic translation initiation factor 4B (eIF4B). Phosphorylation of S6 leads to an increase in protein synthesis and cell proliferation.
基因ID	6199
基因名	RPS6KB2
Swiss	Q9UBS0 ( <a href="https://www.uniprot.org/uniprotkb/Q9UBS0/entry">https://www.uniprot.org/uniprotkb/Q9UBS0/entry</a> )
别名	RPS6KB2,KLS,P70-beta,P70-beta-1,P70-beta-2,S6K-beta2,S6K2,SRK,STK14B,p70(S6K)-beta,p70 S6 Kinase beta Antibody,70 kDa ribosomal protein S6 kinase 2,S6 kinase-related kinase,Serine/threonine-protein kinase 14B,p70 ribosomal S6 kinase beta

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)