

— ABLYBIO, Help Your Research



PKG2 Antibody

货号: **AYP4834**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	84kDa/87kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Rat brain
细胞定位	Apical cell membrane,Lipid-anchor
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human PKG2.
------	--

靶点信息

研究背景	This gene encodes a protein that belongs to the serine/threonine protein kinase family of proteins. The encoded protein plays a role in the regulation of fluid balance in the intestine. A similar protein in mouse is thought to regulate differentiation and proliferation of cells in the colon. Alternate splicing results in multiple transcript variants.
基因ID	5593
基因名	PRKG2
Swiss	Q13237 (https://www.uniprot.org/uniprotkb/Q13237/entry)
别名	PRKGR2,cGK2,cGKII,PRKG2,PKG2 Antibody,cGMP-dependent protein kinase II,PKG2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)