

— ABLYBIO, Help Your Research



TIE2 Antibody

货号: **AYP4842**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	109kDa/121kDa/125kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HUVEC
细胞定位	Cell junction,Cell membrane,Cytoplasm,Secreted,Single-pass type I membrane protein,cytoskeleton,focal adhesion
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human TIE2.
------	--

靶点信息

研究背景	This gene encodes a receptor that belongs to the protein tyrosine kinase Tie2 family. The encoded protein possesses a unique extracellular region that contains two immunoglobulin-like domains, three epidermal growth factor (EGF)-like domains and three fibronectin type III repeats. The ligand angiopoietin-1 binds to this receptor and mediates a signaling pathway that functions in embryonic vascular development. Mutations in this gene are associated with inherited venous malformations of the skin and mucous membranes. Alternative splicing results in multiple transcript variants. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known.
基因ID	7010
基因名	TEK
Swiss	Q02763 (https://www.uniprot.org/uniprotkb/Q02763/entry)
别名	TEK,CD202B,GLC3E,TIE-2,TIE2,VMCM,VMCM1,TIE2 Antibody,Endothelial tyrosine kinase,Tunica interna endothelial cell kinase,Tyrosine kinase with Ig and EGF homology domains-2,Tyrosine-protein kinase receptor TEK,Tyrosine-protein kinase receptor TIE-2,p140 TEK

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)