

ERK1/2 Antibody

货号: **AYP4925**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	36kDa/41kDa/38kDa/40kDa/43kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3
细胞定位	caveola,cytoplasm,cytoskeleton,cytosol,early endosome,endoplasmic reticulum lumen,extracellular region,focal adhesion,Golgi apparatus,late endosome,microtubule organizing center,mitochondrion,mitotic spindle,nucleoplasm,nucleus,plasma membrane
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ERK1/2.
------	--

靶点信息

研究背景	<p>This gene encodes a member of the MAP kinase family. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. The activation of this kinase requires its phosphorylation by upstream kinases. Upon activation, this kinase translocates to the nucleus of the stimulated cells, where it phosphorylates nuclear targets. One study also suggests that this protein acts as a transcriptional repressor independent of its kinase activity. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. Two alternatively spliced transcript variants encoding the same protein, but differing in the UTRs, have been reported for this gene. [provided by RefSeq, Jan 2014]</p>
------	--

基因ID	5594,5595
基因名	MAPK1,MAPK3
Swiss	P28482/P27361
别名	MAPK1/MAPK3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn