

— ABLYBIO, Help Your Research



# ROCK2 Antibody

货号: **AYP4927**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	160kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,THP-1, Mouse lung,Mouse brain,Rat brain
细胞定位	Cell membrane,Cytoplasm,Nucleus,Peripheral membrane protein,centrosome,cytoskeleton,microtubule or ganizing center
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human ROCK2.
------	---

## 靶点信息

研究背景	The protein encoded by this gene is a serine/threonine kinase that regulates cytokinesis, smooth muscle contraction, the formation of actin stress fibers and focal adhesions, and the activation of the c-fos serum response element. This protein, which is an isozyme of ROCK1 is a target for the small GTPase Rho.
基因ID	9475
基因名	ROCK2
Swiss	O75116 ( <a href="https://www.uniprot.org/uniprotkb/O75116/entry">https://www.uniprot.org/uniprotkb/O75116/entry</a> )
别名	ROCK2,ROCK-II,ROCK2 Antibody,Rho kinase 2,Rho-associated,coiled-coil-containing protein kinase 2,coiled-coil-containing protein kinase II,p164 ROCK-2,KIAA0619

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)