

— ABLYBIO, Help Your Research



# HDAC3 Antibody

货号: **AYP4957**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF:</b> 1:50 - 1:200
理论分子量	48kDa/49kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,MCF7,NIH/3T3,C6
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human HDAC3.
------	---

## 靶点信息

研究背景	Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene belongs to the histone deacetylase/acuc/alpha family. It has histone deacetylase activity and represses transcription when tethered to a promoter. It may participate in the regulation of transcription through its binding with the zinc-finger transcription factor YY1. This protein can also down-regulate p53 function and thus modulate cell growth and apoptosis. This gene is regarded as a potential tumor suppressor gene.
基因ID	8841
基因名	HDAC3
Swiss	O15379 ( <a href="https://www.uniprot.org/uniprotkb/O15379/entry">https://www.uniprot.org/uniprotkb/O15379/entry</a> )
别名	HD3,RPD3,RPD3-2,HDAC3,HDAC3 Antibody,Protein deacetylase HDAC3,Protein deacylase HDAC3,SMAP45

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))