

— ABLYBIO, Help Your Research



Keratin 8 Antibody

货号: **AYP4958**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IHC ELISA |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 |
| 理论分子量 | 53kDa/56kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HeLa,Jurkat,Rat liver |
| 细胞定位 | Cytoplasm,Nucleus,Nucleus matrix,nucleoplasm |
| 纯化 | Affinity purification |

抗原信息

| | |
|------|---|
| 抗原信息 | Synthesized peptide derived from Human Keratin 8. |
|------|---|

靶点信息

| | |
|-------|--|
| 研究背景 | This gene is a member of the type II keratin family clustered on the long arm of chromosome 12. Type I and type II keratins heteropolymerize to form intermediate-sized filaments in the cytoplasm of epithelial cells. The product of this gene typically dimerizes with keratin 18 to form an intermediate filament in simple single-layered epithelial cells. This protein plays a role in maintaining cellular structural integrity and also functions in signal transduction and cellular differentiation. Mutations in this gene cause cryptogenic cirrhosis. Alternatively spliced transcript variants have been found for this gene. |
| 基因ID | 3856 |
| 基因名 | KRT8 |
| Swiss | P05787 (https://www.uniprot.org/uniprotkb/P05787/entry) |
| 别名 | KRT8,CARD2,CK-8,CK8,CYK8,K2C8,K8,KO,keratin 8,Keratin 8 Antibody,Cytokeratin-8,Keratin-8,Type-II keratin Kb8 |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)