

— ABLYBIO, Help Your Research



# PIAS1 Antibody

货号: **AYP4968**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	36kDa/71kDa/72kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus,Nucleus speckle,PML body
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human PIAS1.
------	-----------------------------------------------

## 靶点信息

研究背景	This gene encodes a member of the protein inhibitor of activated STAT (PIAS) family. PIAS proteins function as SUMO E3 ligases and play important roles in many cellular processes by mediating the sumoylation of target proteins. This protein plays a central role as a transcriptional coregulator of numerous cellular pathways including the STAT1 and nuclear factor kappaB pathways. Alternate splicing results in multiple transcript variants.
基因ID	8554
基因名	PIAS1
Swiss	O75925 ( <a href="https://www.uniprot.org/uniprotkb/O75925/entry">https://www.uniprot.org/uniprotkb/O75925/entry</a> )
别名	PIAS1,DDXBP1,GBP,GU/RH-II,ZMIZ3,PIAS1 Antibody,DEAD/H box-binding protein 1,E3 SUMO-protein transferase PIAS1,Gu-binding protein,Protein inhibitor of activated STAT protein 1,RNA helicase II-binding protein

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)