

— ABLYBIO, Help Your Research



IF2B3 Antibody

货号: **AYP4971**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	21kDa/63kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,BxPC-3,SGC-7901,NIH/3T3,293T,HepG2
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human IF2B3.
------	---

靶点信息

研究背景	The protein encoded by this gene is primarily found in the nucleolus, where it can bind to the 5' UTR of the insulin-like growth factor II leader 3 mRNA and may repress translation of insulin-like growth factor II during late development. The encoded protein contains several KH domains, which are important in RNA binding and are known to be involved in RNA synthesis and metabolism. A pseudogene exists on chromosome 7, and there are putative pseudogenes on other chromosomes.
基因ID	10643
基因名	IGF2BP3
Swiss	O00425 (https://www.uniprot.org/uniprotkb/O00425/entry)
别名	IGF2BP3,CT98,IMP-3,IMP3,KOC,KOC1,VICKZ3,IF2B3 Antibody,IGF-II mRNA-binding protein 3,KH domain-containing protein overexpressed in cancer,VICKZ family member 3,IF2B3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)