

PFDN1 Antibody

货号: **AYP4977**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	14kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,PC-12,Mouse brain
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human PFDN1.
------	---

靶点信息

研究背景	This gene encodes a member of the prefoldin beta subunit family. The encoded protein is one of six subunits of prefoldin, a molecular chaperone complex that binds and stabilizes newly synthesized polypeptides, thereby allowing them to fold correctly. The complex, consisting of two alpha and four beta subunits, forms a double beta barrel assembly with six protruding coiled-coils.
基因ID	5201
基因名	PFDN1
Swiss	O60925

别名	PFDN1;PDF;PFD1
----	----------------

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn