

— ABLYBIO, Help Your Research



TUSC2 Antibody

货号: **AYP4983**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	12kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	mitochondrion
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human TUSC2.
------	---

靶点信息

研究背景	This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available.
------	--

基因ID	11334
基因名	TUSC2
Swiss	O75896 (https://www.uniprot.org/uniprotkb/O75896/entry)
别名	TUSC2,C3orf11,FUS1,PAP,PDAP2,TUSC2 Antibody,Fusion 1 protein,PDGFA-associated protein 2,LGCC

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)