

— ABLYBIO, Help Your Research



GAS41 Antibody

货号: **AYP4995**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	26kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human GAS41.
------	---

靶点信息

研究背景	The protein encoded by this gene is found in the nucleoli. It has high sequence homology to human MLLT 1, and yeast and human MLLT3 proteins. Both MLLT1 and MLLT3 proteins belong to a class of transcription factors, indicating that the encoded protein might also represent a transcription factor. This protein is thought to be required for RNA transcription. This gene has been shown to be amplified in tumors. Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	8089
基因名	YEATS4
Swiss	O95619 (https://www.uniprot.org/uniprotkb/O95619/entry)
别名	YEATS4,4930573H17Rik,B230215M10Rik,GAS41,NUBI-1,YAF9,GAS41 Antibody,Glioma-amplified sequence 41,NuMA-binding protein 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)