

— ABLYBIO, Help Your Research



TNFA Antibody

货号: **AYP5004**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200
理论分子量	25kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PMA, LPS and Brefeldin A treated THP-1 cells
细胞定位	Cell membrane,Membrane,Secreted,Single-pass type II membrane protein,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human TNFA.
------	--

靶点信息

研究背景	This gene encodes a multifunctional proinflammatory cytokine that belongs to the tumor necrosis factor (TNF) superfamily. This cytokine is mainly secreted by macrophages. It can bind to, and thus functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/TNFR. This cytokine is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation, apoptosis, lipid metabolism, and coagulation. This cytokine has been implicated in a variety of diseases, including autoimmune diseases, insulin resistance, and cancer. Knockout studies in mice also suggested the neuroprotective function of this cytokine.
基因ID	7124
基因名	TNF
Swiss	P01375 (https://www.uniprot.org/uniprotkb/P01375/entry)
别名	TNF,DIF,TNF-alpha,TNFA,TNFSF2,TNLG1F,TNF alpha,TNFA Antibody,Cachectin,Tumor necrosis factor ligand superfamily member 2,N-terminal fragment

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)