

— ABLYBIO, Help Your Research



TACD2 Antibody

货号: **AYP5017**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,A-431,Mouse kidney
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human TACD2.
------	---

靶点信息

研究背景	This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like cornea I dystrophy.
基因ID	4070
基因名	TACSTD2
Swiss	P09758 (https://www.uniprot.org/uniprotkb/P09758/entry)
别名	TACSTD2,EGP-1,EGP1,GA733-1,GA7331,GP50,M1S1,TROP2,TACD2 Antibody,Cell surface glycoprotein Trop-2,Membrane component chromosome 1 surface marker 1,Pancreatic carcinoma marker protein GA733-1,TACD2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)