

— ABLYBIO, Help Your Research



# ATF7 Antibody

货号: **AYP5048**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	13kDa/33kDa/49kDa/50kDa/51kDa/52kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,SW620,HepG2
细胞定位	Cytoplasm,Nucleus,nucleoplasm
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human ATF7.
------	--

## 靶点信息

研究背景	Stress-responsive chromatin regulator that plays a role in various biological processes including innate immunological memory, adipocyte differentiation or telomerase regulation. In absence of stress, contributes to the formation of heterochromatin and heterochromatin-like structure by recruiting histone H3K9 tri- and di-methyltransferases thus silencing the transcription of target genes such as STAT1 in adipocytes, or genes involved in innate immunity in macrophages and adipocytes (By similarity. Stress induces ATF7 phosphorylation that disrupts interactions with histone methyltransferase and enhances the association with coactivators containing histone acetyltransferase and/or histone demethylase, leading to disruption of the heterochromatin-like structure and subsequently transcriptional activation (By similarity. In response to TNF-alpha, which is induced by various stresses, phosphorylated ATF7 and telomerase are released from telomeres leading to telomere shortening. Plays also a role in maintaining epithelial regenerative capacity and protecting against cell death during intestinal epithelial damage and repair (By similarity.
基因ID	11016
基因名	ATF7
Swiss	P17544 ( <a href="https://www.uniprot.org/uniprotkb/P17544/entry">https://www.uniprot.org/uniprotkb/P17544/entry</a> )
别名	ATF7,ATFA,ATF7 Antibody,Activating transcription factor 7,Transcription factor ATF-A

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))