

— ABLYBIO, Help Your Research



# EGR2 Antibody

货号: **AYP5056**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	44kDa/50kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse kidney
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human EGR2.
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a transcription factor with three tandem C2H2-type zinc fingers. Defects in this gene are associated with Charcot-Marie-Tooth disease type 1D (CMT1D), Charcot-Marie-Tooth disease type 4E (CMT4E), and with Dejerine-Sottas syndrome (DSS). Multiple transcript variants encoding two different isoforms have been found for this gene.
基因ID	1959
基因名	EGR2
Swiss	P11161 ( <a href="https://www.uniprot.org/uniprotkb/P11161/entry">https://www.uniprot.org/uniprotkb/P11161/entry</a> )
别名	EGR2,AT591,CMT1D,CMT4E,KROX20,EGR2 Antibody,E3 SUMO-protein transferase ERG2,Early growth response protein 2,Zinc finger protein Krox-20

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))