

TNF Receptor II Antibody

货号: **AYP5065**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IF ELISA |
| 推荐浓度 | WB: 1:500 - 1:2000 IF: 1:50 - 1:200 |
| 理论分子量 | 28kDa/48kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | SW620,Mouse liver,Mouse spleen,Rat liver |
| 细胞定位 | Cell membrane,Secreted,Single-pass type I membrane protein |
| 纯化 | Affinity purification |

抗原信息

| | |
|------|---|
| 抗原信息 | Synthesized peptide derived from Human TNF Receptor II. |
|------|---|

靶点信息

| | |
|------|---|
| 研究背景 | The protein encoded by this gene is a member of the TNF-receptor superfamily. This protein and TNF-receptor 1 form a heterocomplex that mediates the recruitment of two anti-apoptotic proteins, c-IAP1 and c-IAP2, which possess E3 ubiquitin ligase activity. The function of IAPs in TNF-receptor signalling is unknown, however, c-IAP1 is thought to potentiate TNF-induced apoptosis by the ubiquitination and degradation of TNF-receptor-associated factor 2, which mediates anti-apoptotic signals. Knockout studies in mice also suggest a role of this protein in protecting neurons from apoptosis by stimulating antioxidative pathways. |
| 基因ID | 7133 |
| 基因名 | TNFRSF1B |

| | |
|--------------|--|
| Swiss | P20333 |
| 别名 | CD120b;TBPII;TNF-R-II;TNF-R75;TNFBR;TNFR1B;TNFR2;TNFR80;p75;p75TNFR;TNF Receptor II;TNFRSF1B |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn