

— ABLYBIO, Help Your Research



GRK2 Antibody

货号: **AYP5078**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	79kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Mouse spleen,Mouse brain
细胞定位	Cell membrane,Cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human GRK2.
------	--

靶点信息

研究背景	The product of this gene phosphorylates the beta-2-adrenergic receptor and appears to mediate agonist-specific desensitization observed at high agonist concentrations. This protein is an ubiquitous cytosolic enzyme that specifically phosphorylates the activated form of the beta-adrenergic and related G-protein-coupled receptors. Abnormal coupling of beta-adrenergic receptor to G protein is involved in the pathogenesis of the failing heart.
基因ID	156
基因名	GRK2
Swiss	P25098 (https://www.uniprot.org/uniprotkb/P25098/entry)
别名	GRK2,ADRBK1,BARK1,BETA-ARK1,GRK2 Antibody,G-protein coupled receptor kinase 2,BARK

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)