

— ABLYBIO, Help Your Research



# MARK3 Antibody

货号: **AYP5083**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	human brain tissue, PC-3 cells
细胞定位	cytoplasm,cytosol,dendrite,extracellular exosome,plasma membrane
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human MARK3.
------	---

## 靶点信息

研究背景	The protein encoded by this gene is activated by phosphorylation and in turn is involved in the phosphorylation of tau proteins MAP2 and MAP4. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]
------	--

基因ID	4140
基因名	MARK3
Swiss	P27448 ( <a href="https://www.uniprot.org/uniprotkb/P27448/entry">https://www.uniprot.org/uniprotkb/P27448/entry</a> )
别名	CTAK1,KP78,PAR1A,Par-1a,MARK3,MARK3 Antibody,C-TAK1,Cdc25C-associated protein kinase 1,ELKL motif kinase 2,Protein kinase STK10,Ser/Thr protein kinase PAR-1,Serine/threonine-protein kinase p78,EMK2

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)