

— ABLYBIO, Help Your Research



# CORO1A Antibody

货号: AYP5098

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	51kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,SW480,Mouse brain,Mouse lung,Rat lung
细胞定位	Cytoplasm,Cytoplasmic vesicle,cell cortex,cytoskeleton,phagosome membrane
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human CORO1A.
------	--

## 靶点信息

研究背景	This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. Alternative splicing results in multiple transcript variants. A related pseudogene has been defined on chromosome 16.
基因ID	11151
基因名	CORO1A
Swiss	P31146 ( <a href="https://www.uniprot.org/uniprotkb/P31146/entry">https://www.uniprot.org/uniprotkb/P31146/entry</a> )
别名	CORO1A,CLABP,CLIPINA,HCORO1,IMD8,TACO,p57,CORO1A Antibody,Coronin-like protein A,Coronin-like protein p57,Tryptophan aspartate-containing coat protein,CORO1

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)