

— ABLYBIO, Help Your Research



RORA Antibody

货号: **AYP5108**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200
理论分子量	53kDa/58kDa/61kDa/63kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human RORA.
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The encoded protein has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Also, it has been shown to aid in the transcriptional regulation of some genes involved in circadian rhythm. Four transcript variants encoding different isoforms have been described for this gene.
基因ID	6095
基因名	RORA
Swiss	P35398 (https://www.uniprot.org/uniprotkb/P35398/entry)
别名	RORA,NR1F1,ROR1,ROR2,ROR3,RZR-ALPHA,RZRA,RORA Antibody,Nuclear receptor RZR-alpha,Nuclear receptor subfamily 1 group F member 1,RAR-related orphan receptor A,Retinoid-related orphan receptor-alpha

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)