

— ABLYBIO, Help Your Research



# GRK6 Antibody

货号: **AYP5120**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	64kDa/65kDa/67kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Jurkat,293T
细胞定位	Lipid-anchor,Membrane
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human GRK6.
------	--

## 靶点信息

研究背景	This gene encodes a member of the guanine nucleotide-binding protein (G protein)-coupled receptor kinase subfamily of the Ser/Thr protein kinase family. The protein phosphorylates the activated forms of G protein-coupled receptors thus initiating their deactivation. Several transcript variants encoding different isoforms have been described for this gene.
基因ID	2870
基因名	GRK6
Swiss	P43250 ( <a href="https://www.uniprot.org/uniprotkb/P43250/entry">https://www.uniprot.org/uniprotkb/P43250/entry</a> )
别名	GRK6,GPRK6,GRK6 Antibody,G protein-coupled receptor kinase GRK6

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn//www.ablybio.cn>)