

— ABLYBIO, Help Your Research



# CLK1 Antibody

货号: **AYP5135**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	16kDa/57kDa/61kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Jurkat
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human CLK1.
------	--

## 靶点信息

研究背景	This gene encodes a member of the CDC2-like (or LAMMER) family of dual specificity protein kinases. In the nucleus, the encoded protein phosphorylates serine/arginine-rich proteins involved in pre-mRNA processing, releasing them into the nucleoplasm. The choice of splice sites during pre-mRNA processing may be regulated by the concentration of transacting factors, including serine/arginine rich proteins. Therefore, the encoded protein may play an indirect role in governing splice site selection. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	1195
基因名	CLK1
Swiss	P49759 ( <a href="https://www.uniprot.org/uniprotkb/P49759/entry">https://www.uniprot.org/uniprotkb/P49759/entry</a> )
别名	CLK1,CLK,CLK/STY,STY,CLK1 Antibody,CDC-like kinase 1

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn//www.ablybio.cn>)