

— ABLYBIO, Help Your Research



# p57 Kip2 Antibody

货号: **AYP5139**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	31kDa/32kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse lung,Rat kidney
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human p57 Kip2.
------	--

## 靶点信息

研究背景	This gene is imprinted, with preferential expression of the maternal allele. The encoded protein is a tight-binding, strong inhibitor of several G1 cyclin/Cdk complexes and a negative regulator of cell proliferation. Mutations in this gene are implicated in sporadic cancers and Beckwith-Wiedemann syndrome, suggesting that this gene is a tumor suppressor candidate. Three transcript variants encoding two different isoforms have been found for this gene.
基因ID	1028
基因名	CDKN1C
Swiss	P49918 ( <a href="https://www.uniprot.org/uniprotkb/P49918/entry">https://www.uniprot.org/uniprotkb/P49918/entry</a> )
别名	CDKN1C,BWCR,BWS,KIP2,WBS,p57,p57Kip2,p57 Kip2 Antibody,Cyclin-dependent kinase inhibitor p57,p57 Kip2

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn//www.ablybio.cn>)