

— ABLYBIO, Help Your Research



GALK1 Antibody

货号: **AYP5145**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	42kDa/45kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HepG2,HeLa,MCF7,Mouse heart,Mouse brain,Mouse lung
细胞定位	cytoplasm,cytosol,extracellular exosome
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human GALK1.
------	---

靶点信息

研究背景	Galactokinase is a major enzyme for the metabolism of galactose and its deficiency causes congenital cataracts during infancy and presenile cataracts in the adult population.
------	--

基因ID	2584
基因名	GALK1
Swiss	P51570 (https://www.uniprot.org/uniprotkb/P51570/entry)
别名	GALK1,GALK,GK1,HEL-S-19,GALK1 Antibody,Galactose kinase

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)