

— ABLYBIO, Help Your Research



MKK6Antibody

货号: **AYP5156**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	31kDa/37kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cytoplasm,Nucleus,cytoskeleton
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human MKK6Antibody.
------	------------------------------------------------------

靶点信息

研究背景	This gene encodes a member of the dual specificity protein kinase family, which functions as a mitogen-activated protein (MAP) kinase kinase. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals. This protein phosphorylates and activates p38 MAP kinase in response to inflammatory cytokines or environmental stress. As an essential component of p38 MAP kinase mediated signal transduction pathway, this gene is involved in many cellular processes such as stress induced cell cycle arrest, transcription activation and apoptosis.
基因ID	5608
基因名	MAP2K6
Swiss	P52564 (https://www.uniprot.org/uniprotkb/P52564/entry)
别名	MAP2K6,MAPKK6,MEK6,MKK6,PRKMK6,SAPKK-3,SAPKK3,MKK6Antibody,MAPK/ERK kinase 6,Stress-activated protein kinase kinase 3,SKK3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)