

— ABLYBIO, Help Your Research



TPD52 Antibody

货号: **AYP5177**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200
理论分子量	16kDa/19kDa/22kDa/24kDa/25kDa/26kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,HeLa,HT-29,A375
细胞定位	cytoplasm,endoplasmic reticulum,perinuclear region of cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human TPD52.
------	---

靶点信息

研究背景	Enables calcium ion binding activity and protein homodimerization activity. Involved in B cell differentiation. Located in endoplasmic reticulum and perinuclear region of cytoplasm.
基因ID	7163
基因名	TPD52
Swiss	P55327 (https://www.uniprot.org/uniprotkb/P55327/entry)
别名	TPD52,D52,N8L,PC-1,PrLZ,hD52,TPD52 Antibody,Protein N8

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)