

— ABLYBIO, Help Your Research



# MtSSB Antibody

货号: **AYP5215**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF:</b> 1:50 - 1:200
理论分子量	17kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,HL-60,293T,Mouse liver
细胞定位	Mitochondrion,Mitochondrion matrix,mitochondrion nucleoid
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human MtSSB.
------	---

## 靶点信息

研究背景	SSBP1 is a housekeeping gene involved in mitochondrial biogenesis (Tiranti et al., 1995 [PubMed 778999 1]). It is also a subunit of a single-stranded DNA (ssDNA)-binding complex involved in the maintenance of genome stability (Huang et al., 2009) [PubMed 19683501].[supplied by OMIM, Feb 2010]
基因ID	6742
基因名	SSBP1
Swiss	Q04837 ( <a href="https://www.uniprot.org/uniprotkb/Q04837/entry">https://www.uniprot.org/uniprotkb/Q04837/entry</a> )
别名	SSBP1,Mt-SSB,SOSS-B1,SSBP,mtSSB,MtSSB Antibody,PWP1-interacting protein 17

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))