

— ABLYBIO, Help Your Research



# CSTF1 Antibody

货号: **AYP5216**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	48kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U2OS,MCF7,SW480,Mouse thymus,Rat thymus
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human CSTF1.
------	---

## 靶点信息

研究背景	This gene encodes one of three subunits which combine to form cleavage stimulation factor (CSTF). CSTF is involved in the polyadenylation and 3'end cleavage of pre-mRNAs. Similar to mammalian G protein bet a subunits, this protein contains transducin-like repeats. Several transcript variants with different 5' UTR, but encoding the same protein, have been found for this gene.
基因ID	1477
基因名	CSTF1
Swiss	Q05048 ( <a href="https://www.uniprot.org/uniprotkb/Q05048/entry">https://www.uniprot.org/uniprotkb/Q05048/entry</a> )
别名	CSTF1,CstF-50,CstFp50,CSTF1 Antibody,CF-1 50 kDa subunit,Cleavage stimulation factor 50 kDa subunit

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)