

— ABLYBIO, Help Your Research



SFRS4 Antibody

货号: **AYP5221**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	56kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat uterus,Rat thymus
细胞定位	Nucleus speckle
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human SFRS4.
------	---

靶点信息

研究背景	This gene encodes a member of the arginine/serine-rich splicing factor family. The encoded protein likely functions in mRNA processing.
------	---

基因ID	6429
基因名	SRSF4
Swiss	Q08170 (https://www.uniprot.org/uniprotkb/Q08170/entry)
别名	SRSF4,SFRS4,SRP75,SFRS4 Antibody,Pre-mRNA-splicing factor SRP75,SRP001LB,Splicing factor,arginine/serine-rich 4

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)