

— ABLYBIO, Help Your Research



KLF11 Antibody

货号: **AYP5239**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse thymus,mouse brain, rat lung
细胞定位	cytosol,focal adhesion,nuclear body,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human KLF11.
------	---

靶点信息

研究背景	
------	--

基因ID	8462
基因名	KLF11
Swiss	O14901 (https://www.uniprot.org/uniprotkb/O14901/entry)
别名	KLF11 Antibody,KLF11,Transforming growth factor-beta-inducible early growth response protein 2,FKLF,TI EG2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)