

— ABLYBIO, Help Your Research



# TCOF1 Antibody

货号: **AYP5247**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF:</b> 1:50 - 1:200
理论分子量	96kDa/144kDa/148kDa/152kDa/155kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,293T
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human TCOF1.
------	---

## 靶点信息

研究背景	This gene encodes a nucleolar protein with a LIS1 homology domain. The protein is involved in ribosomal DNA gene transcription through its interaction with upstream binding factor (UBF). Mutations in this gene have been associated with Treacher Collins syndrome, a disorder which includes abnormal craniofacial development. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	6949
基因名	TCOF1
Swiss	Q13428 ( <a href="https://www.uniprot.org/uniprotkb/Q13428/entry">https://www.uniprot.org/uniprotkb/Q13428/entry</a> )
别名	TCOF1,MFD1,TCS,TCS1,treacle,TCOF1 Antibody,Treacher Collins syndrome protein

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)