

— ABLYBIO, Help Your Research



# PTK6/BRK Antibody

货号: **AYP5261**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	14kDa/51kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,BT-474,MCF7,A-549,NIH/3T3,Mouse liver,Mouse kidney,Mouse testis,Rat kidney
细胞定位	Cell projection,Cytoplasm,Membrane,Nucleus,ruffle
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human PTK6/BRK.
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a cytoplasmic nonreceptor protein kinase which may function as an intracellular signal transducer in epithelial tissues. Overexpression of this gene in mammary epithelial cells leads to sensitization of the cells to epidermal growth factor and results in a partially transformed phenotype. Expression of this gene has been detected at low levels in some breast tumors but not in normal breast tissue. The encoded protein has been shown to undergo autophosphorylation. Alternative splicing results in multiple transcript variants.
基因ID	5753
基因名	PTK6
Swiss	Q13882 ( <a href="https://www.uniprot.org/uniprotkb/Q13882/entry">https://www.uniprot.org/uniprotkb/Q13882/entry</a> )
别名	PTK6, BRK, PTK6/BRK Antibody, Breast tumor kinase, Tyrosine-protein kinase BRK, PTK6/BRK

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)