

— ABLYBIO, Help Your Research



# ANKRD1 Antibody

货号: **AYP5286**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	36kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human ANKRD1.
------	--

## 靶点信息

研究背景	The protein encoded by this gene is localized to the nucleus of endothelial cells and is induced by IL-1 and TNF-alpha stimulation. Studies in rat cardiomyocytes suggest that this gene functions as a transcription factor. Interactions between this protein and the sarcomeric proteins myopalladin and titin suggest that it may also be involved in the myofibrillar stretch-sensor system.
基因ID	27063
基因名	ANKRD1
Swiss	Q15327 ( <a href="https://www.uniprot.org/uniprotkb/Q15327/entry">https://www.uniprot.org/uniprotkb/Q15327/entry</a> )
别名	ANKRD1,ALRP,C-193,CARP,CVARP,MCARP,bA320F15.2,ANKRD1 Antibody,Cardiac ankyrin repeat protein,Cytokine-inducible gene C-193 protein,Cytokine-inducible nuclear protein,C193,HA1A2

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)