

— ABLYBIO, Help Your Research



AKAP14 Antibody

货号: **AYP5351**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	23kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,DU145,Mouse testis,Rat testis
细胞定位	axoneme
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human AKAP14.
------	--

靶点信息

研究背景	The A-kinase anchor proteins (AKAPs) are a group of structurally diverse proteins, which have the common function of binding to the regulatory subunit of protein kinase A (PKA) and confining the holoenzyme to discrete locations within the cell. This gene encodes a member of the AKAP family. The protein anchors PKA in ciliary axonemes and, in this way, may play a role in regulating ciliary beat frequency. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]
基因ID	158798
基因名	AKAP14
Swiss	Q86UN6 (https://www.uniprot.org/uniprotkb/Q86UN6/entry)
别名	AKAP28,PRKA14,AKAP14 Antibody,AKAP14,A-kinase anchor protein 28 kDa,Protein kinase A-anchoring protein 14

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)