

— ABLYBIO, Help Your Research



CAMK5 Antibody

货号: **AYP5371**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200
理论分子量	51kDa/54kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse eye,Mouse spleen,Rat eye
细胞定位	glutamatergic synapse,plasma membrane
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human CAMK5.
------	---

靶点信息

研究背景	Does not appear to have detectable kinase activity.
------	---

基因ID	79012
基因名	CAMKV
Swiss	Q8NCB2 (https://www.uniprot.org/uniprotkb/Q8NCB2/entry)
别名	CAMKV,1G5,VACAMKL,CAMK5,CAMK5 Antibody

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)